

Anti-SARS-CoV-2 NSP14 antibody (380-stop) (STJ140192)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Goat
Reactivity	SARS-CoV-2

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	3 mg/mL
Conjugation	Unconjugated
Purification	Epitope affinity purified
Dilution Range	WB 1:500-1:2000
Formulation	PBS, 20% Glycerol and 0.05% Sodium Azide.
Isotype	IgG
Storage Instruction	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.

TARGET INFORMATION

Immunogen Region	380-stop
Immunogen Sequence	MVCLFWNCNVD RYPANSIVC VCLFWNCNVD RYPANSIVCR CLFWNCNVD RYPANSIVCRF LFWNCNVD RYPANSIVCRFD FWNCNVD RYPANSIVCRFDT WNCNVD RYPANSIVCRFDTR NCNVD RYPANSIVCRFDTRV CNVD RYPANSIVCRFDTRVL
Specificity	In lysates of transfected cells with the plasmid containing the sequence used, detects the fusion protein by Western blot.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Anti SARS-CoV-2 NSP14 Goat Polyclonal
antibody western blot STJ140192