

## Goat Anti-Human IgM (Full length) antibody {Unconjugated} (STJ140199)

### GENERAL INFORMATION

Product Type	Secondary antibodies
Applications	WB
Host / Source	Goat
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	2 mg/mL
Conjugation	Unconjugated
Purification	Epitope affinity purified
Dilution Range	WB 1:1000-1:5000
Formulation	PBS, 20% Glycerol and 0.05% Sodium Azide.
Isotype	IgG
Storage Instruction	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.

### TARGET INFORMATION

Immunogen Region	Full length
Specificity	Using serum samples detects the light and heavy chains by Western blot.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---