

## Goat Anti-Rat IgG (Full length) antibody {DyLight®488} (STJ140235)

### GENERAL INFORMATION

Product Type	Secondary antibodies
Applications	WB/IF
Host / Source	Goat
Reactivity	Rat

### PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	2 mg/mL
Conjugation	DyLight®488
Purification	Epitope affinity purified
Dilution Range	WB 1:2000-1:10000 IF 1:100-1:1000
Formulation	PBS, 20% Glycerol, 0.05% ProClin® and 0.05% Sodium Azide.
Isotype	IgG
Storage Instruction	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.

### TARGET INFORMATION

Immunogen Region	Full length
Specificity	Using serum samples detects rat IgG by Western blot. Reduced cross-reactivity to other species.

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Anti IgG Goat Polyclonal antibody  
 immunofluorescence mouse rat STJ140235