

Anti-CANX antibody (550 aa C-Term) {DyLight®488} (STJ140256)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF/IHC-F
Host / Source	Goat
Reactivity	Human/Rat/Mouse/Canine/Monkey

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	2 mg/mL
Conjugation	DyLight®488
Purification	Epitope affinity purified
Dilution Range	WB 1:500-1:5000 IF 1:50-1:500 IHC-F 1:200-1:1000
Formulation	PBS, 20% Glycerol, 0.05% ProClin® and 0.05% Sodium Azide.
Isotype	IgG
Storage Instruction	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.

TARGET INFORMATION


Gene ID	821
Gene Symbol	CANX
UniProt ID	CALX_HUMAN
Immunogen Region	550 aa C-Term
Immunogen Sequence	KSDAEEDGGTVSQEEEDRKP SDAEEDGGTVSQEEEDRKP
Specificity	Detects a band of 90 kDa by Western blot whole cell lysates.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Anti CANX Goat Polyclonal antibody western blot
hela mcf7 STJ140256



Anti CANX Goat Polyclonal antibody
immunofluorescence STJ140256