

Anti-GAPDH antibody (240 aa C-Term) {DyLight®405} (STJ140264)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	WB/IF/IHC-F
Host / Source	Goat
Reactivity	Human/Rat/Mouse/Zebrafish/Canine/Monkey
PRODUCT PROPERTIES	
Clonality	Polyclonal
Concentration	2 mg/mL
Conjugation	DyLight®405
Purification	Epitope affinity purified
Dilution Range	WB 1:500-1:5000 IHC-F 1:200-1:1000 IF 1:50-1:250
Formulation	PBS, 20% Glycerol, 0.05% ProClin® and 0.05% Sodium Azide.
Isotype	IgG
Storage Instruction	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.
TARGET INFORMATION	
Gene ID	2597
Gene Symbol	GAPDH
UniProt ID	G3P_HUMAN
Immunogen Region	240 aa C-Term
Immunogen Sequence	SVVDLTCRLEKPAKYDDIKK VVDLTCRLEKPAKYDDIKKV VDLTCRLEKPAKYDDIKKV DLTCRLEKPAKYDDIKKVVK LTCRLEKPAKYDDIKKVVKQ
Specificity	Detects a band of 37 kDa by Western blot.
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.