

Anti-GAPDH antibody (240 aa C-Term) {DyLight®550} (STJ140266)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF/IHC-F
Host / Source	Goat
Reactivity	Human/Rat/Mouse/Zebrafish/Canine/Monkey

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	2 mg/mL
Conjugation	DyLight®550
Purification	Epitope affinity purified
Dilution Range	WB 1:500-1:5000 IHC-F 1:200-1:1000 IF 1:50-1:250
Formulation	PBS, 20% Glycerol, 0.05% ProClin® and 0.05% Sodium Azide.
Isotype	IgG
Storage Instruction	For continuous use, store at 2-8 C for one-two days. For extended storage, store in-20 C freezer. Working dilution samples should be discarded if not used within 12 hours.

TARGET INFORMATION

Gene ID	2597
Gene Symbol	GAPDH
UniProt ID	G3P_HUMAN
Immunogen Region	240 aa C-Term
Immunogen Sequence	SVVDLTCRLEKPAKYDDIKK VVDLTCRLEKPAKYDDIKKV VDLTCRLEKPAKYDDIKKVV DLTCRLEKPAKYDDIKKVVK LTCRLEKPAKYDDIKKVVKQ
Specificity	Detects a band of 37 kDa by Western blot.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Anti GAPDH Goat Polyclonal antibody western blot cos7 jurkat STJ140266