

## Anti-mGrape antibody (STJ140281)

### GENERAL INFORMATION

|               |                       |
|---------------|-----------------------|
| Product Type  | Primary antibodies    |
| Applications  | WB/IF/IHC-P/IHC-F/IEM |
| Host / Source | Goat                  |
| Reactivity    | RFP                   |

### PRODUCT PROPERTIES


|                     |                                                                                                                                                                              |
|---------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Clonality           | Polyclonal                                                                                                                                                                   |
| Concentration       | 3 mg/mL                                                                                                                                                                      |
| Conjugation         | Unconjugated                                                                                                                                                                 |
| Purification        | Epitope affinity purified                                                                                                                                                    |
| Dilution Range      | WB 1:500-1:5000<br>IHC-F 1:50-1:500<br>IHC-P 1:50-1:500<br>IF 1:50-1:500<br>IEM 1:50-1:500                                                                                   |
| Formulation         | PBS, 20% Glycerol and 0.05% Sodium Azide.                                                                                                                                    |
| Isotype             | IgG                                                                                                                                                                          |
| Storage Instruction | For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours. |

### TARGET INFORMATION


|             |                                                                                                                                                            |
|-------------|------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Specificity | In 293HEK cells transfected with cds plasmid detects a band of 25.4 kDa by Western blot. This antibody does not recognize GFP (green fluorescent protein). |
|-------------|------------------------------------------------------------------------------------------------------------------------------------------------------------|

### ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Anti mGrape Goat Polyclonal  
immunofluorescence STJ140281



Anti mGrape Goat Polyclonal antibody western  
blot hek293 STJ140281



Anti mGrape Goat Polyclonal antibody  
immunofluorescence STJ140281