

Anti-CD9 antibody (115-195 aa) (STJ140294)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Goat
Reactivity	Human/Rat/Mouse/Canine/Monkey

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	3 mg/mL
Conjugation	Unconjugated
Purification	Epitope affinity purified
Dilution Range	WB:1:1, 000-1:2, 000
Formulation	PBS, 20% Glycerol and 0.05% Sodium Azide.
Isotype	IgG
Storage Instruction	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.

TARGET INFORMATION

Immunogen Region	115-195 aa
Immunogen Sequence	HKDEVIKEVQEFYKDTYNKL KDEVIKEVQEFYKDTYNKLK DEVIKEVQEFYKDTYNKLKT EVIKEVQEFYKDTYNKLKTK VIKEVQEFYKDTYNKLKTKD
Specificity	Reacts with CD9 detected by Western blot using the following human (Jurkat, MCF7, HeLa, HEK-293) , mouse (NIH3T3) whole cell lysates.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Anti CD9 Goat Polyclonal antibody western blot
hek293 hela STJ140294