

Anti-NusA antibody (STJ140298)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Goat
Reactivity	E.coli

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	3 mg/mL
Conjugation	Unconjugated
Purification	Epitope affinity purified
Dilution Range	WB:1:1, 000-1:5, 000
Formulation	PBS, 20% Glycerol and 0.05% Sodium Azide.
Isotype	IgG
Storage Instruction	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.

TARGET INFORMATION

Immunogen Sequence	MGSSHHHHHHSSMNKEILAV GSSHHHHHHSSMNKEILAVV SSHHHHHHHSSMNKEILAVVE SHHHHHHHSSMNKEILAVVEA HHHHHHHSSMNKEILAVVEAV HHHHHSSMNKEILAVVEAVS HHHHSSMNKEILAVVEAVSN HHHSSMNKEILAVVEAVSNE HHHSSMNKEILAVVEAVSNEK HSSMNKEILAVVEAVSNEKA SSMNKEILAVVEAVSNEKAL SMNKEILAVVEAVSNEKAL
Specificity	Detects NusA protein by Western blot.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Anti NusA Goat Polyclonal antibody western blot
STJ140298