

Anti-GAPDH antibody (C-Term) {DyLight®680} (STJ140310)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF/IHC-F
Host / Source	Goat
Reactivity	Human/Rat/Mouse/Zebrafish/Canine/Monkey

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	2.5 mg/mL
Conjugation	DyLight®680
Purification	Epitope affinity purified
Dilution Range	WB:1:500-1:5,000 IF:1:50-1:250 IHC-F:1:200-1:1,000
Formulation	PBS, 20% Glycerol and 0.05% Sodium Azide.
Isotype	IgG
Storage Instruction	For continuous use, store at 2-8 C for one-two days. For extended storage, store in-20 C freezer. Working dilution samples should be discarded if not used within 12 hours.

TARGET INFORMATION

Immunogen Region	C-Term
Immunogen Sequence	SVVDLTCRLEKPAKYDDIKK VVDLTCRLEKPAKYDDIKKV VDLTCRLEKPAKYDDIKKVV DLTCRLEKPAKYDDIKKVVK LTCRLEKPAKYDDIKKVVKQ
Specificity	Detects a band of 37 kDa by Western blot.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Anti GAPDH Goat Polyclonal antibody western blot hela jurkat STJ140310