

Anti-NPY3R antibody (C-Term) (STJ140317)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Goat
Reactivity	Human/Rat/Mouse/Canine/Monkey

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	3 mg/mL
Conjugation	Unconjugated
Purification	Epitope affinity purified
Dilution Range	WB:1:500-1:2, 000
Formulation	PBS, 20% Glycerol and 0.05% Sodium Azide.
Isotype	IgG
Storage Instruction	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.

TARGET INFORMATION

Immunogen Region	C-Term
Immunogen Sequence	GAKFKTSAQHALTSVSRGSS AKFKTSAQHALTSVSRGSSL KFKTSAQHALTSVSRGSSLK
Specificity	Detects endogenous levels of NPY3R by Western blot. This Ab is specific for NPY3R and does not cross-react with NPY1R, NPY2R, NPY4R and NPY5R.

ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



Anti NPY3R Goat Polyclonal antibody western blot jurkat STJ140317