

## Anti-LacZ antibody (STJ140318)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Goat
Reactivity	E.coli

### PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	3 mg/mL
Conjugation	Unconjugated
Purification	Epitope affinity purified
Dilution Range	WB:1:500-1:5, 000
Formulation	PBS, 20% Glycerol and 0.05% Sodium Azide.
Isotype	IgG
Storage Instruction	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.

### TARGET INFORMATION

Immunogen Sequence	MTMITDSLAVVLQRRDWENP TMITDSLAVVLQRRDWENPG MITDSLAVVLQRRDWENPGV ITDSLAVVLQRRDWENPGVT TDSLAVVLQRRDWENPGVTQ DSLAVVLQRRDWENPGVTQL SLAVVLQRRDWENPGVTQLN LAVVLQRRDWENPGVTQLNR AVVLQRRDWENPGVTQLNRL VVLQRRDWENPGVTQLNRLA VLQRRDWENPGVTQLNRLAA LQRRDWENPGVTQLNRLAA
Specificity	In lysates of transfected/transformed cells with the plasmid containing the sequence used, detects the fusion protein by Western blot.

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Anti LacZ Goat Polyclonal antibody western blot  
hela STJ140318