

## Mouse IgG1 Isotype Control [NA] {PE} (STJ16100001)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.1 mg/mL
Conjugation	PE
Formulation	PBS with 0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

### TARGET INFORMATION

Gene Symbol	<a href="#">NAUniprot ID=IGG1_HUMAN"</a>
-------------	--

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---