

Mouse IgG2a Isotype Control [NA] {FITC} (STJ16100002)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	50 tests
Conjugation	FITC
Formulation	PBS with 0.1% sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---