

## Anti-CD8 antibody [FK18] (STJ16100010)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IF
Host / Source	Mouse
Reactivity	Human/Rhesus

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Dilution Range	IF 1:10
Formulation	Serum free medium containing 1% FCS
Isotype	IgG3
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---