

Anti-CD57 antibody [6-13-19-1] (STJ16100011)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IF/IHC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Dilution Range	IF 1:5 IHC 1:5
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgM
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---