

Anti-CD11c antibody [FK24] (STJ16100012)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IF/IHC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Dilution Range	IF 1:5 IHC 1:5
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	3687
Gene Symbol	ITGAX
UniProt ID	ITAX_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---