

Anti-CD11a antibody [87-6F9] (STJ16100022)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/FUNC/WB/IHC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution Range	Flow cytometry (0, 5-1, 0 µg/million cells in 0, 1 ml). Functional studies (0, 02-2, 0 µg/ml without azide). Immunoblotting (1-2 µg/ml). Immunohistology (1-2 µg/ml for 30 min at RT; an appropriate antigen retrieval method for staining of formalin-fix
Formulation	PBS with 0.02% Sodium Azide.
Isotype	IgG2bk
Storage Instruction	Store for up to 1 year at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	3683
Gene Symbol	ITGAL
UniProt ID	ITAL_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---