

Anti-Beta-2 microglobulin antibody [B2M-01] (STJ16100027)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IP/WB/IHC-P/ICC/ELISA
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Dilution Range	FC 0 5-1 0 µg/million cells in 0 1 ml.-WB 2-4 µg/ml-ELISA to be determined-IHC paraffin Recommended
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

