

Anti-CD37 antibody [IPO-24] (STJ16100030)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IHC-F/IF
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	100 µg/mL
Conjugation	Unconjugated
Dilution Range	FC 0 5-1 0 µg/million cells in 0 1 ml.-IF 0 5-1 0 µg/ml.-IHC 1-2 µg/ml for 30 min at RT
Formulation	PBS with 0.02% sodium azide
Isotype	IgG2b
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	951
Gene Symbol	CD37
UniProt ID	CD37_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---