

Anti-CD29 antibody [MEM-101A] (STJ16100040)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IP/WB
Host / Source	Mouse
Reactivity	Dog/Pig/Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Dilution Range	FC 4 µg/ml-IP to be determined
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

