

## Anti-CD28 antibody [204.12] (STJ16100044)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IHC-F/IF
Host / Source	Mouse
Reactivity	Human/Mouse/Pig/Sheep/Cow

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	100 $\mu$ g/mL
Conjugation	Unconjugated
Dilution Range	FC 0 5-1 0 $\mu$ g/million cells in 0 1 ml-IF 0 5-1 0 $\mu$ g/ml-IHC frozen 1-2 $\mu$ g/ml for 30-60 min at RT
Formulation	PBS with 0.02% sodium azide
Isotype	IgG2b
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

### TARGET INFORMATION

Gene ID	<a href="#">940</a>
Gene Symbol	<a href="#">CD28</a>
UniProt ID	<a href="#">CD28_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---