

## Anti-HLA-DR+DP antibody [HL-40] (STJ16100047)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IP/IHC-P
Host / Source	Mouse
Reactivity	Human/Pig

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Dilution Range	FC 1 $\mu$ g/ml-IP-IHC paraffin to be determined
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

