

## Anti-HLA-DQ1+DQ3 antibody [HL-37] (STJ16100048)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | FC/IP/WB           |
| Host / Source | Mouse              |
| Reactivity    | Pig/Human          |

### PRODUCT PROPERTIES

|                     |  |
|---------------------|--|
| Clonality           | Monoclonal   |
| Concentration       | 1 mg/mL  |
| Conjugation         | Unconjugated   |
| Dilution Range      | FC 5 µg/ml-IP-WB to be determined                                |
| Formulation         | Phosphate buffered saline (PBS) solution with 15 mM sodium azide |
| Isotype             | IgG3   |
| Storage Instruction | Store at 2-8°C for up to 1-year, upon receipt.                   |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

