

## Anti-CD4 antibody [MEM-115] (STJ16100050)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IP
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Dilution Range	FC 3 µg/ml-IP
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

### TARGET INFORMATION

Gene ID	<a href="#">920</a>
Gene Symbol	<a href="#">CD4</a>
UniProt ID	<a href="#">CD4_HUMAN</a>

### ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**

