

## Anti-CD14 antibody [MEM-18] (STJ16100056)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IP/WB/ELISA
Host / Source	Mouse
Reactivity	Human/Non-Human Primates

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Dilution Range	FC 4 $\mu$ g/ml-IP-WB Non-reducing conditions.-ELISA The antibody MEM-18 has been tested as the detection antibody
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

### TARGET INFORMATION

Gene ID	<a href="#">929</a>
Gene Symbol	<a href="#">CD14</a>
UniProt ID	<a href="#">CD14_HUMAN</a>

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

