

Anti-CD44 antibody [MEM-85] {PE} (STJ16100064)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	PE
Dilution Range	FC 1:50
Formulation	PBS with sucrose, BSA and 0.1% sodium azide
Isotype	IgG2b
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	960
Gene Symbol	CD44
UniProt ID	CD44_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---