

Anti-PECAM-1 antibody [MEC7.46] (STJ16100083)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	IHC-F/IHC-P/ELISA/FC/IF/IP
Host / Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	100 Mu g/mL
Conjugation	Unconjugated
Dilution Range	IHC-F 1:50 IF 1:50 FC 1:50
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---