

Anti-CD55 antibody [MEM-118] (STJ16100088)

GENERAL INFORMATION

| | |
|---------------|--------------------------|
| Product Type | Primary antibodies |
| Applications | FC/IP/IHC-P/WB |
| Host / Source | Mouse |
| Reactivity | Human/Non-Human Primates |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Clonality | Monoclonal |
| Concentration | 1 mg/mL |
| Conjugation | Unconjugated |
| Formulation | Tris buffered saline (TBS) solution with 15 mM sodium azide |
| Isotype | IgM |
| Storage Instruction | Store at 2-8°C for up to 1-year, upon receipt. |

TARGET INFORMATION

| | |
|-------------|---------------------------|
| Gene ID | 1604 |
| Gene Symbol | CD55 |
| UniProt ID | DAF_HUMAN |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|