

Anti-CD34 antibody [QBEND-10] (STJ16100095)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/IHC-P/WB/FC/IP/IF/EM
Host / Source	Mouse
Reactivity	Human/Rhesus/Cynomolgus Monkey

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1.0 mg/mL
Conjugation	Unconjugated
Dilution Range	FC 1:100-1:200 IHC 1:500-1:1000
Formulation	PBS with 0.09% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	947
Gene Symbol	CD34
UniProt ID	CD34_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---