

Anti-CD34 antibody [QBEND-10] {Biotin} (STJ16100096)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Human/Rhesus/Cynomolgus Monkey

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.1mg/mL
Conjugation	Biotin
Formulation	PBS with 0.09% sodium azide and 1% BSA
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	947
Gene Symbol	CD34
UniProt ID	CD34_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---