

Anti-CD47 antibody [BRIC126] (STJ16100098)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/FC/IHC-F/IP
Host / Source	Mouse
Reactivity	Human/Pig/Sheep/Dog/Bovine

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Dilution Range	FC 1:10
Formulation	TRIS buffered glycine
Isotype	IgG2b
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	961
Gene Symbol	CD47
UniProt ID	CD47_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---