

Anti-CD51 antibody [23C6] (STJ16100100)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IF/IP/FC/IHC-F
Host / Source	Mouse
Reactivity	Human/Chicken

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1.0 mg/mL
Conjugation	Unconjugated
Dilution Range	FC 1:40-1:80
Formulation	PBS with 0.09% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---