

Anti-CD16 antibody [3G8] {PE} (STJ16100112)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Human/Rhesus

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	PE
Formulation	PBS with 4mg/ml BSA and 0.1% sodium azide
Isotype	IgG1k
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---