

Anti-CD20 antibody [HI47] {PE} (STJ16100115)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Human/Rhesus

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	PE
Formulation	PBS with 4mg/ml BSA, sucrose and 0.1% sodium azide
Isotype	IgG3
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	931
Gene Symbol	MS4A1
UniProt ID	CD20_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---