

Anti-HLA-Class I antibody [W6/32] (STJ16100127)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IP/WB/IHC-F/ICC/ELISA
Host / Source	Mouse
Reactivity	Cow/Non-Human Primates/Human/Cat

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---