

Anti-Nucleophosmin antibody [FC-61991] (STJ16100133)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/ICC/IF/ELISA/IP
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Dilution Range	ELISA 1-2 µg/mL ICC 2-10 µg/mL IF 2-10 µg/mL
Formulation	PBS with 0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	4869
Gene Symbol	NPM1
UniProt ID	NPM_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---