

Anti-S. epidermis antibody [17-5] (STJ16100140)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA
Host / Source	Mouse
Reactivity	Bacteria

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	200 Mu g/mL
Conjugation	Unconjugated
Dilution Range	ELISA 1:10
Formulation	RPMI with 0.02% sodium azide
Isotype	IgG3
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---