

Anti-Lipid A antibody [26-5] (STJ16100142)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/AGG
Host / Source	Mouse
Reactivity	Bacteria

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	100 Mu g/mL
Conjugation	Unconjugated
Formulation	PBS with 0.02% sodium azide and 0.1% BSA
Isotype	IgG2b
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---