

Anti-CD100 antibody [133.1G6] (STJ16100144)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IHC-F
Host / Source	Mouse
Reactivity	Human/Cow/Baboon

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	100 Mu g/mL
Conjugation	Unconjugated
Formulation	PBS with 0.02% sodium azide
Isotype	IgM
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	10507
Gene Symbol	SEMA4D
UniProt ID	SEM4D_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---