

## Anti-B3 Integrin/CD61 antibody [BV4] (STJ16100154)

### GENERAL INFORMATION

|               |                         |
|---------------|-------------------------|
| Product Type  | Primary antibodies      |
| Applications  | IHC-F/ELISA/IP/IHC-P/WB |
| Host / Source | Mouse                   |
| Reactivity    | Human/Bovine            |

### PRODUCT PROPERTIES

|                     |  |
|---------------------|--|
| Clonality           | Monoclonal                                     |
| Concentration       | 100 Mu g/mL                                    |
| Conjugation         | Unconjugated                                   |
| Dilution Range      | FC 1:50<br>ELISA 1:50<br>WB 1:50<br>IP 1:50    |
| Formulation         | PBS with 0.1% BSA and 0.02% sodium azide       |
| Isotype             | IgG1   |
| Storage Instruction | Store at 2-8°C for up to 1-year, upon receipt. |

### ADDITIONAL INFORMATION

|      |   |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|