

## Anti-Myeloid Specific Antigen antibody [BM-2] (STJ16100158)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IHC
Host / Source	Mouse
Reactivity	Human/Macaque Monkey

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution Range	Flow cytometry (0.5-1.0 µg/million cells in 0.1 ml). Immunohistology (1-2 µg/ml for 30 min at RT; staining of formalin-fixed tissues better after boiling tissue sections in 10mM Citrate Buffer, pH6, 0, for 10-20 min followed by cooling at RT for 20 m
Formulation	PBS with 0.02% Sodium Azide.
Isotype	IgG1
Storage Instruction	Store for up to 1 year at 2-8°C upon receipt.

### ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**

