

## Anti-Blood Group A (CD173) antibody [HE-10] (STJ16100160)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	AGG/FC/IF/IHC
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution Range	Agglutination. Flow cytometry (0, 5–1, 0 µg/million cells in 0, 1 ml). Immunofluorescence (0, 5–1 µg/ml). Immunohistology (formalin-fixed: 1-2 µg/ml for 30 min. at RT; requires boiling tissue sections in 10mM citrate buffer, pH 6, 0, for 10-20 min fo
Formulation	PBS with 0.02% Sodium Azide.
Isotype	IgMk
Storage Instruction	Store for up to 1 year at 2-8°C upon receipt.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---