

## Anti-blood group antigen A antibody [HE-195] (STJ16100164)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FUNC
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Formulation	Tris buffered saline (TBS) solution with 15 mM sodium azide
Isotype	IgM
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

### ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**

Clone	Agglutination of tested blood groups				
	A <sub>1</sub>	B	A-B	A <sub>2</sub> B	0
HE-193 (anti-A)	+	-	+	+	-
HE-195 (anti-A)	+	-	+	+	-
HEB-23 (anti-B)	-	+	+	+	-
HE-24 (anti-A:B)	+	-	+	-	-