

Anti-NCAM antibody [MEM-188] (STJ16100190)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IP/IHC-P/ICC
Host / Source	Mouse
Reactivity	Non-Human Primates/Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

