

Anti-Mannose Receptor antibody [15-2] (STJ16100193)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | IHC-F/FC/WB |
| Host / Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Clonality | Monoclonal |
| Concentration | 100 Mu g/mL |
| Conjugation | Unconjugated |
| Dilution Range | WB 1:50 |
| Formulation | PBS with 0.1% BSA and 0.02% sodium azide |
| Isotype | IgG1 |
| Storage Instruction | Store at 2-8°C for up to 1-year, upon receipt. |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|