

Anti-HLA-E antibody [MEM-E/06] (STJ16100204)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IP/IHC-P
Host / Source	Mouse
Reactivity	Non-Human Primates/Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	3133
Gene Symbol	HLA-E
UniProt ID	HLAE_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

