

Anti-STAT 1 antibody [SM1] (STJ16100213)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IP/WB/IHC-P
Host / Source	Mouse
Reactivity	Mouse/Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG2b
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**

