

Anti-Nectin-2 antibody [502-57] (STJ16100216)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/FC/ELISA/IP/WB
Host / Source	Rat
Reactivity	Mouse/Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	100 Mu g/mL
Conjugation	Unconjugated
Dilution Range	FC 1:50 IHC 1:50 ELISA 1:50 WB 1:50 IP 1:50
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	5819
Gene Symbol	NECTIN2
UniProt ID	NECT2_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---