

Anti-TLR9 antibody [5G5] (STJ16100217)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/FC/IF/IHC-P/WB
Host / Source	Mouse
Reactivity	Human/Dog/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	100 Mu g/mL
Conjugation	Unconjugated
Dilution Range	FC 1:50 IHC 1:50 WB 1:50 IF 1:50
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	54106
Gene Symbol	TLR9
UniProt ID	TLR9_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---